

Quest 88 limited

Aston Street, SHIFNAL, Shropshire TF11 8DW
 tel: +44 (0)1952 463050 fax: +44 (0)1952 463077

USER INSTRUCTIONS Upi Multi-Stander conversion to hip stander.

These are the instructions for converting an Upi Multi-stander into a hip stander. Please read them carefully before attempting to convert the stander. Please contact Quest 88 if you are unsure as to how to proceed with the procedure. Do not attempt this procedure when a user is standing in the frame.

The standard Upi has hip and thoracic supports. If a user has relatively good upper trunk control or the user's therapist wishes to work on upper trunk exercises, the multi stander can be converted to a simple hip stander.

When this is done, the lower hip support of the stander is removed and the main height adjustment tubes shortened so that the thoracic support can be set at a low height, suitable for supporting the user at the hips.

To convert the stander pull out the circle clip pins which control height adjustment on the frame (see main user instructions). Lift out the whole support section.

Undo the knobs on the sides of the hip support and slide the hip support off the height adjustment tubes.

Another set of circle clip pins fit into the height adjustment tubes below the thoracic support. Pull these out and remove the extension tubes.

Fit the thoracic support back onto the Upi base to use the Upi as a hip stander.

Reverse this process to convert the Upi back to a full vertical stander.

See the main Upi instructions for operation and set up of the frame.

